

UT Southwestern Medical Center at Dallas
Southwestern Center for Minimally Invasive Surgery

Trainee Name: _____

FLS Training Program
Task 1: Peg Transfer

Proficiency Score: 48 seconds with no pegs dropped out of the field of view

Proficiency Training Protocol: Achieve proficiency score on 2 consecutive repetitions and on 10 additional (nonconsecutive is acceptable) reps for reinforcement ("**2 + 10**")
If proficiency is not reached by 80 repetitions, proceed to next task

Trainee Scores

Date	Repetition No.	Time	No. Pegs Not Transferred
	1		
	2		
	3		
	4		
	5		
	6		
	7		
	8		
	9		
	10		
	11		
	12		
	13		
	14		
	15		
	16		
	17		
	18		
	19		
	20		
	21		
	22		
	23		
	24		
	25		
	26		
	27		
	28		
	29		
	30		
	31		
	32		
	33		
	34		
	35		
	36		
	37		
	38		
	39		
	40		

Date	Repetition No.	Time	No. Pegs Not Transferred
	41		
	42		
	43		
	44		
	45		
	46		
	47		
	48		
	49		
	50		
	51		
	52		
	53		
	54		
	55		
	56		
	57		
	58		
	59		
	60		
	61		
	62		
	63		
	64		
	65		
	66		
	67		
	68		
	69		
	70		
	71		
	72		
	73		
	74		
	75		
	76		
	77		
	78		
	79		
	80		

UT Southwestern Medical Center at Dallas
Southwestern Center for Minimally Invasive Surgery

Trainee Name: _____

FLS Training Program
Task 2: Pattern Cut

Proficiency Score: 98 seconds with all cuts within 2mm of the line (either side)

Proficiency Training Protocol: Achieve proficiency score on 2 consecutive repetitions only ("**2 Consecutive**")
If proficiency is not reached by 80 repetitions, proceed to next task

Trainee Scores

Date	Repetition No.	Time	Outside of 2mm Accuracy (Yes or No)
	1		
	2		
	3		
	4		
	5		
	6		
	7		
	8		
	9		
	10		
	11		
	12		
	13		
	14		
	15		
	16		
	17		
	18		
	19		
	20		
	21		
	22		
	23		
	24		
	25		
	26		
	27		
	28		
	29		
	30		
	31		
	32		
	33		
	34		
	35		
	36		
	37		
	38		
	39		
	40		

Date	Repetition No.	Time	Outside of 2mm Accuracy (Yes or No)
	41		
	42		
	43		
	44		
	45		
	46		
	47		
	48		
	49		
	50		
	51		
	52		
	53		
	54		
	55		
	56		
	57		
	58		
	59		
	60		
	61		
	62		
	63		
	64		
	65		
	66		
	67		
	68		
	69		
	70		
	71		
	72		
	73		
	74		
	75		
	76		
	77		
	78		
	79		
	80		

UT Southwestern Medical Center at Dallas
Southwestern Center for Minimally Invasive Surgery

Trainee Name: _____

FLS Training Program
Task 3: Endoloop

Proficiency Score: 53 seconds with up to 1mm accuracy error; no knot security errors (slippage) allowed

Proficiency Training Protocol: Achieve proficiency score on 2 consecutive repetitions only ("**2 Consecutive**")
If proficiency is not reached by 80 repetitions, proceed to next task

Trainee Scores

Date	Repetition No.	Time	Outside of 1mm Accuracy (Yes or No)	Knot Security Error (Yes or No)
	1			
	2			
	3			
	4			
	5			
	6			
	7			
	8			
	9			
	10			
	11			
	12			
	13			
	14			
	15			
	16			
	17			
	18			
	19			
	20			
	21			
	22			
	23			
	24			
	25			
	26			
	27			
	28			
	29			
	30			
	31			
	32			
	33			
	34			
	35			
	36			
	37			
	38			
	39			
	40			

Date	Repetition No.	Time	Outside of 1mm Accuracy (Yes or No)	Knot Security Error (Yes or No)
	41			
	42			
	43			
	44			
	45			
	46			
	47			
	48			
	49			
	50			
	51			
	52			
	53			
	54			
	55			
	56			
	57			
	58			
	59			
	60			
	61			
	62			
	63			
	64			
	65			
	66			
	67			
	68			
	69			
	70			
	71			
	72			
	73			
	74			
	75			
	76			
	77			
	78			
	79			
	80			

UT Southwestern Medical Center at Dallas
Southwestern Center for Minimally Invasive Surgery

Trainee Name: _____

FLS Training Program
Task 4: Extracorporeal Suture

Proficiency Score: 136 seconds with up to 1mm accuracy error; no knot security errors (slippage) allowed

Proficiency Training Protocol: Achieve proficiency score on 2 consecutive repetitions only ("**2 Consecutive**")
If proficiency is not reached by 80 repetitions, proceed to next task

Trainee Scores

Date	Repetition No.	Time	Outside of 1mm Accuracy (Yes or No)	Knot Security Error (Yes or No)
	1			
	2			
	3			
	4			
	5			
	6			
	7			
	8			
	9			
	10			
	11			
	12			
	13			
	14			
	15			
	16			
	17			
	18			
	19			
	20			
	21			
	22			
	23			
	24			
	25			
	26			
	27			
	28			
	29			
	30			
	31			
	32			
	33			
	34			
	35			
	36			
	37			
	38			
	39			
	40			

Date	Repetition No.	Time	Outside of 1mm Accuracy (Yes or No)	Knot Security Error (Yes or No)
	41			
	42			
	43			
	44			
	45			
	46			
	47			
	48			
	49			
	50			
	51			
	52			
	53			
	54			
	55			
	56			
	57			
	58			
	59			
	60			
	61			
	62			
	63			
	64			
	65			
	66			
	67			
	68			
	69			
	70			
	71			
	72			
	73			
	74			
	75			
	76			
	77			
	78			
	79			
	80			

UT Southwestern Medical Center at Dallas
Southwestern Center for Minimally Invasive Surgery

Trainee Name: _____

FLS Training Program
Task 5: Intracorporeal Suture

Proficiency Score: 112 seconds with up to 1mm accuracy error; no knot security errors (slippage) allowed

Proficiency Training Protocol: Achieve proficiency score on 2 consecutive repetitions and on 10 additional (nonconsecutive is acceptable) reps for reinforcement ("**2 + 10**")

If proficiency is not reached by 80 repetitions, proceed to next task

Trainee Scores

Date	Repetition No.	Time	Outside of 1mm Accuracy (Yes or No)	Knot Security Error (Yes or No)
	1			
	2			
	3			
	4			
	5			
	6			
	7			
	8			
	9			
	10			
	11			
	12			
	13			
	14			
	15			
	16			
	17			
	18			
	19			
	20			
	21			
	22			
	23			
	24			
	25			
	26			
	27			
	28			
	29			
	30			
	31			
	32			
	33			
	34			
	35			
	36			
	37			
	38			
	39			
	40			

Date	Repetition No.	Time	Outside of 1mm Accuracy (Yes or No)	Knot Security Error (Yes or No)
	41			
	42			
	43			
	44			
	45			
	46			
	47			
	48			
	49			
	50			
	51			
	52			
	53			
	54			
	55			
	56			
	57			
	58			
	59			
	60			
	61			
	62			
	63			
	64			
	65			
	66			
	67			
	68			
	69			
	70			
	71			
	72			
	73			
	74			
	75			
	76			
	77			
	78			
	79			
	80			

